



MILK DOCTOR TETRA

- ▶ Unnecessary extraction the samples
- ▶ Specifically match to MRL level
- ▶ Result in 5 minutes
- ▶ All materials in one kit
- ▶ Simple operation and user friendly



MADE
IN
JAPAN

50test

SIMPLE!

EASY!

FASTNESS!

USER FRIENDLY

Unnecessary to
extract samples

NO INSTRUMENTS

Unnecessary
incubator to heat samples

HIGH SENSITIVITY

Tetracycline 20ppb
Oxytetracycline 20ppb

EYES INTERPRETATION

Easy to get result


Typical Tetracycline Detection Limits

- Tetracycline	20ppb
- Oxytetracycline	20ppb
- Chlortetracycline	100ppb
- Doxycycline	70ppb



Operating Procedure


Collect milk
samples



Add milk
samples



Dip test paper and
wait at room
temperature



	TC	CTC	OTC	DC	
0ppb	50	100	50	50	0ppb

AR Brown Co., Ltd Specialty Chemical dept.

Agri Food safety FD-B team

Dai-ni Marutaka Bldg., 7-13-8, Ginza, Chuo-ku, Tokyo, Japan

TEL : 81-3-3545-5724 FAX : 81-3-3543-8865

URL : <http://www.arbrown.com/> MAIL: a-r@arbrown.com